

**FORCEPS** بینس  
Ophthalmic Surgical Instruments

**www.binasaz.com**

8/0



**Tying Forceps**  
Straight, with micro platforms  
8 cm long  
B# 7835  
scale : 1


10/0



**Tying Forceps**  
Straight, with micro platforms  
10 cm long  
B# 7850C  
scale : 1



10/0



**Tying Forceps**  
Straight, with micro platforms  
10 cm long (curved)  
B# 13168  
scale : 1



**Jewelers Forceps**  
B# 3050  
scale : 1



**Castroviejo Forceps**  
 Oblique teeth 0.4 mm, and  
 tying platform  
 B# 7348  
 scale : 1



**Castroviejo Forceps**  
 Oblique teeth 0.12 mm, and  
 tying platform  
 B# 3045  
 scale : 1



**Fixtation Forceps 1\*2 teeth**  
 B# 3225  
 scale : 1



**Akahoshi Phaco PreChopper**  
 For cracking, dividing and rotating  
 the nucleus without sculpting  
 B# As-130  
 B# As-140 (curved)  
 scale : 1





**Barraquer-Troutman forceps**  
"Colibri" type, oblique teeth 0.12 mm  
with tying platforms  
B# 7862



**Mc Pherson Forceps**  
Ergonomic, angled  
platform 11 mm  
very thin. 12 cm long  
B# 13242



**Mc Pherson Forceps**  
Long smooth jaw  
platform 11 mm  
B# 17161





**Capsulorhexis Forceps**  
Ergonomic, curved  
cystotome type  
11.5 cm long  
B# 18259  
scale : 1

8/0



**Castroviejo Needle Holder**  
Delicate, Stright jaws,without lock  
gentle spring action  
B# 8148  
scale : 1

10/0



**Barraquer-De Laage Needle Holder**  
Curved jaws 1 mm, without lock  
13 cm long  
B# 7251  
scale : 1